

TASMAN/E • TASMAN/N • TASMAN/EN

Hydrotherapy Systems



(ordering code 01.501.080 • 01.501.081 • 01.501.078) Four-Cell Bath

Electrotherapeutic treatments in the water produce beneficial effects on various systems (nervous, circulatory) of the human body. The Four-Cell Bath allows treatment to be concentrated on a specific area of the body (lower or upper limb, one side of the body). According to the indications, the therapist can set the polarity of the electrode located in each cell and choose between galvanic, diadynamic and faradic current. All parameters are entered via a user-friendly electronic control panel. Alternating treatments use water to induce thermal stimuli on the limbs submerged in upper and lower cells with the purpose of stimulating reactions in the circulatory system. The duration of the stimulation, the water temperature and the number of cycles of cold/warm water can be set by therapist.

TASMAN – The Four-Cell Bath is equipped with hand shower for easy cleaning the tub after the treatment. The bathtub has new design and ergonomic construction.

TASMAN - The Four-Cell Bath allows three kinds of treatment to be performed in the water environment:

- **Tasman/E** (electrotherapy treatment)
- **Tasman/N** (alternate treatment)
- **Tasman/EN** (combination of the above)



SPOLKA Z O.O.
TECHNOMEX

MAIN FEATURES

Advantages of TASMAN:

- Spoutless water filling system
- Thermomixer facilitating regulation of water temperature during filling
- Clear and intelligible electronic control panel
- Modern design and ergonomic construction

STANDARD EQUIPMENT FOR TASMAN/EN

- Hand shower
- Filling system
- Electronic control panel for electrotherapy and alternate water treatment settings
- Stainless steel electrodes for electrotherapy treatment
- Valve for adjusting water pressure for alternate treatment
- Thermomixer facilitating regulation of water temperature during filling
- Universal chair BETA with the height regulation: 54-79 cm [1] - (19.220.001)



STANDARD AVAILABLE COLOURS

	5005 RAL Signal Blue
	5012 RAL Light Blue
	6027 RAL Light Green
	6018 RAL Yellow Green
	6002 RAL Leaf Green
	10217 RAL White
	9001 RAL Cream
	3017 RAL Rose
	2004 RAL Pure Orange

Specify the colour inside and outside:

- Standard colour is blue (5012 RAL) inside.
- Standard colour is white (10217 RAL) outside.

SPECIFICATIONS

Model of bathtub	Type	Working capacity of chamber (l)		External dimensions (mm) (LxWxH)
		upper	lower	
01.501.080	TASMAN/E	15	25	1240 x 1050 x 1000
01.501.081	TASMAN/N	15	25	1240 x 1050 x 1000
01.501.078	TASMAN/EN	15	25	1240 x 1050 x 1000

