

KOLUMB • KOLUMB BASIC

(ordering code 01.501.052 • 01.501.053) Lower extremities and lumbar spine whirlpool tub

Lower extremities and lumbar spine whirlpool **KOLUMB** is very popular model of whirlpool for rehabilitation centers, fitness clubs and wellness centers. KOLUMB tub is used for whirlpool treatments of lower extremities and lumbar spine, pearl bath and chromotherapy.

Ability to perform treatment both for lower extremities and lumbar spine, low costs and low water capacity (as compared to full body underwater massage tubs) make this product very useful not only for clinics but also for wellness centers.

KOLUMB tub is dedicated for whirlpool massage and consists of a system of wall-mounted 14 jets, divided into 3 independent sections with outflow water stream control. System of jets function control enables free choice of an area that is intended to be treated and also combined function of jets (7 possibilities) for example ankles and lumbar spine. Moreover due to 2 levels of filling (half or full tub) it enables saving water. During the treatment, the patient sits on a seat inside the tub.

Optionally **KOLUMB** can be equipped with coloured touch control panel, which allows for programming:

- Total time of treatment
- Automatic turn off after finish of treatment
- Free programs - user can save 100 programs
- AutoTest - automatic test of safety working fuses



CHROMOTHERAPY SYSTEM
Emission 7 basic colors and mixed colors

SPOLKA Z O.O.
TECHNOMEX

MAIN FEATURES

The biggest advantages of KOLUMB:

- spoutless water filling system
- pearl bath system (55 jets: 11x5 jets)
- chromotherapy
- disinfection system
- AutoTest - automatic test of safety working of fuses - model KOLUMB
- modern design and ergonomic construction

KOLUMB is equipped with modern available systems:

- chemical disinfection - electronic controlled
- automatic filling of tub
- coloured, touch control panel - model KOLUMB

We offer KOLUMB in two versions:

- **KOLUMB BASIC** (electronic control panel)
- **KOLUMB** (coloured, touch control panel)



Electronic control panel - KOLUMB BASIC



Coloured, touch control panel - KOLUMB



Simple and readable software- KOLUMB

STANDARD AVAILABLE COLOURS

	5005 RAL Signal Blue
	5012 RAL Light Blue
	6027 RAL Light Green
	6018 RAL Yellow Green
	6002 RAL Leaf Green
	10217 RAL White
	9001 RAL Cream
	3017 RAL Rose
	2004 RAL Pure Orange

Specify the colour inside and outside:

- Standard colour is blue (5012 RAL) inside.
- Standard colour is white (10217 RAL) outside.



CE 0197



Optionally **KOLUMB** can be equipped with chromotherapy system - 20 multicoloured diodes (7 main colours and derivated colours).

Independent digital panel lets for:

- control of time of lighting of colours, individual programming of cycles of alternative lighting



[1]



SPECIFICATIONS

ELEMENTS OF EQUIPMENT

Type	KOLUMB	KOLUMB BASIC	Elements of equipment	KOLUMB*	KOLUMB BASIC*
Model:	01.501.052	01.501.053	Hand shower	+	+
Working capacity (I) Level I:	80	80	Filling system	automatic	automatic
Working capacity (I) Level II:	210	210	Setting of desired water level	automatic	automatic
External Dimensions (mm)(LxWxH):	1540 x 870 x 1070		Control panel	coloured, touch	electronic
			AutoTest - automatic test of safety working fuses	+	-
			Pump "dry-working" protection	+	+
			Automatic release of water Electronic controlled	option	option
			Timer	+	+
			Automatic switch off after time of treatment	+	+
			Pearl bath system	+	option
			Chromotherapy	+	option
			Disinfection system	+	option
			Anal jets	option	option

STANDARD EQUIPMENT

- Stair 53L x 43W x 39H cm [1] - (19.037.001)



*Optionally the tub can be equipped with more jets, for example for whirlpool massage or pearl bath system.