

(ordering code 08.501.001) Spine traction table

PLATINIUM is the ideal solution of traction treatment, for patients that suffer any lumbar pathology (e.g. discopathy, protrusion and hernia) and cervical pathology (e.g. neck pain, hernia, cervical sprain, etc).

This new traction bed is designed to allow treatment in 3 positions; flexion/extension, rotation and lateral flexion. The possibility to move the patient in three different ways allows the physiotherapist to carry on the treatment with a wide range of positions of the column.

This special bed offer the most refined performance for traction of mobilization/manipulation and new treatment. Its new design and features offer a wide range of treatment for spinal and cervical traction.



Watch the VIDEO!



SPOLKA Z O.O.
TECHNOMEX

MAIN FEATURES

In a very modern and smart design the machine is offering a lot interactive features:

- Performing the traction in static, harmonic, intermittent changing all parameters and storing in the memory.
- Possibility of creating and saving user defined programs.
- Patient database: store patient data with treatment history.
- 7" color touch screen for the easiest and extremely simple use of the machine.
- Smooth adjustment of the rotation and lateral bending of the bed with use of gas springs.
- Adjustable headrest with hole for the face for comfort during the treatment.
- Extremely precise control of the motors, to perform a very comfortable traction.
- Performing the traction treatment in the supine or prone position, in the symmetrical and asymmetrical manner.

TECHNICAL DATA

Power supply:	160 W
Power input:	110-240 V / 50-60 Hz
Weight:	145 Kg
Dimensions (L x W):	227 x 73 cm
Electronical adjustable Height:	58-88 cm
Safety Class:	I B
3 Parts table with:	
- Section for head with adjustable position:	-44° and +36°
- Section for chest	
- Section for pelvic and legs with adjustable position:	
• Headrest:	-44° and +36°
• Flexion/Extension:	-14° and +12°
• Adduction/Abduction:	-20° and +20°
• Internal/External rotat.:	-12° and +12°
- Kinds of traction:	constant and harmonious

- Maximum strenght of the cervical traction:	18 Kg.
- Maximum strenght of lumbar traction:	90 Kg.
- Adjustable following parameters:	type of traction, power, speed, escalation, time of treatment, active phase and breack time
- Possibility to program:	programmable traction phase sequence in a single operation
- Display:	7" color touch screen display
- Maximum load:	150 Kg.

STANDARD EQUIPMENT

- 2 pcs. - Belts [1]: (08.022.003)
 • one to stabilize the lumbar spine
 • one to stabilize the toracic spine
 1 pc. - Stool [2] (08.501.003)



Equipment compliant with
Directive 93/42/EEC integrated by
Directive 2007/47/EC

