


AMETYST

Treatment Tables

(ordering code 13.501.005) Quadriceps Bench

Treatment table/ quadriceps bench that allows exercising and stretching leg muscles. The exercises may be personalized according to patient needs. Load / weight on the lower limbs and angle of training may be adjustable. Painted steel structure, chromed steel leg supports, padded covering, adjustable back rest.



TECHNICAL DATA

Length:	1700 mm
Width:	1000 mm
Height:	760 mm
Height of back rest:	1610 mm
Back rest adjustment:	-30° - +90°

STANDARD EQUIPMENT

- 2 - stabilization belts for leg
- dumb-bells:
 - 2 - 0,5kg
 - 2 - 1,0kg
 - 2 - 1,25 kg



STANDARD AVAILABLE COLOURS

ecological leather



OPTIONAL COLOURS

natural leather

