


# AZURYT 2

Treatment Tables

(ordering code 14.501.009) Treatment Table



This rehabilitation table is designed for use in orthopedic, neurologic, rheumatic diseases treatment as well as in rehabilitation. Places of use are: hospitals, clinics, health centers and rehabilitation centers. The horizontal / vertical position of the table plate is electrically adjusted with the use of remote control.



## TECHNICAL DATA

|                                                    |                                                                        |
|----------------------------------------------------|------------------------------------------------------------------------|
| Power supply:                                      | 230 V - 50 Hz                                                          |
| Power of actuator:                                 | 6000 N                                                                 |
| Power of actuator II (tilting process):            | 8000 N                                                                 |
| As a rule the upholstery consists of two elements: | footrest and main plate                                                |
| Dimensions of the table (L x P x H):               | 200 x 65 x 58-96 cm                                                    |
| Weight of the table (without patient):             | 105 Kg                                                                 |
| Angle adjustment:                                  | - of the main plate +87°<br>- table height adjustment range 58 - 96 cm |
| Maximum patient weight:                            | 150 Kg.                                                                |

## STANDARD AVAILABLE COLOURS

ecological leather



## STANDARD EQUIPMENT

|                                                                     |              |
|---------------------------------------------------------------------|--------------|
| 3 pcs. - Belts for stabilization [1]                                | (13.022.092) |
| 1 pc. - Fixing rods for stabilization belts (head, trunk, head) [2] | (14.501.017) |

## ATTENTION!

Avoid surcharge of the table plate. Maximum surcharge (in elevation and lowering) is 130 Kg, static surcharge is no more than 150 Kg evenly distributed on the whole surface of the table.

## OPTIONAL EQUIPMENT

|                                                       |              |
|-------------------------------------------------------|--------------|
| • Table with adjustable height [3]                    | (14.501.025) |
| • Rechargeable battery                                | (14.501.024) |
| • Central wheel system for mobility (with olders) [4] | (14.501.016) |

## OPTIONAL COLOURS

natural leather

