

TF-204

Treatment Tables

(ordering code 19.022.002) Wooden Table

This **wooden table** is designed for use in physiotherapy, made of high quality wood, doesn't include metal parts.
Places of use are: hospitals, clinics, health centers and rehabilitation centers.



TECHNICAL DATA

Dimensions:	1870mm x 660mm
Maximum weight of patient:	150kg
Weight of table:	30kg
Height:	700mm
Headrest angle-adjustment:	from 0° up to 35°

OPTIONAL EQUIPMENT

• Paper roll holder	19.022.009
---------------------	------------