

MG Wave (evo)

Electro Therapy Equipment



(ordering code 7W00) Magnetic field therapy

MG WAVE is a portable unit for therapy with low-frequency. Pulsating Electromagnetic Fields (ELF-PEMF)

PATIENTS DATABASE

allows storing:

- patients data
- pathology
- description of treatment
- Treatment customisation
- Evaluation of pain level at the beginning and at the end of treatment

Our patient database provides immediate access to the customized treatment allowing it to be easily and quickly repeated.

PROGRAMS AND SEQUENCES

Huge memory to store new programs and sequences of combined programs. this allows creating customized treatments tailored to the patient's specific needs.



MAIN FEATURES

- High output power
- Simple parameter setting
- Pre-set programs for the most common pathologies
- Easy coil positioning
- Flat coils or circular coil for cervical or lumbar treatments
- Ample graphic display showing all output parameters
- Patients database
- Possibility to store modified programs and protocol sequences

STANDARD EQUIPMENT

- 2 pcs. - Ø15 cm flat coils [1] (A01.00083)
 1 pc. - Red cable [2] (A01.00077)
 1 pc. - Black cable [3] (A01.00078)
 2 pcs. - 100 cm elastic straps [4] (F7)
 1 pc. - Magnet [5] (MGN)



Equipment compliant with
 Directive 93/42/EEC integrated by
 Directive 2007/47/EC



OPTIONAL EQUIPMENT

- Ø 30 cm limb coil [6] (A01.00064)
- Ø 60 cm lumbar coil [7] (A01.00063)
- Aluminium therapy bed [8] (A01.00006)
- Wooden therapy bed [9] (A01.00040)
- Red Ø 7.5 cm coil [10] (A01.00071_R)
- Black Ø 7.5 cm coil [11] (A01.00071_N)
- Limb/lumbar coil cable [12] (A01.00065)
- Vest [13] (GIL)
- Evoline bag (A02.00241)
- Evoline trolley [14] (A01.00056)
- Tray for Evoline trolley [15] (A01.00047)



TECHNICAL DATA

Number of output: 2
Field Pilot Signal: rectangular
Pulse Duration: up to 5 ms
Frequency: 1 ÷ 750 Hz
Max. B induction: 160 Gauss (Peak) for Ø 15 cm and Ø 7.5 cm coils 65 Gauss (Peak) for Ø 30 cm coil 50 Gauss (Peak) for Ø 60 cm coil
Timer: 0÷240 min., end-of-treatment acoustic signal and automatic interruption of operation
Power Supply: 115 - 230 Vac, 50 - 60 Hz
Input: 55 VA
Dimensions: - Central unit: 27 x 32.8 x 23 cm (LxWxH) - Flat coils: Ø 15 cm each
Weight: - Central unit: 3.8 Kg - Flat coils: 750 gr. each
Insulation class: Class I - BF type

